

---

# Peritoneal Dialysis

---

## PD Solutions Designed to optimize patient prescription

---



### **Dianeal®**

For acute or chronic renal failure patients

PD4: Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

PD101: Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L

Available in 0.5, 1.5, 2.5 and 4.25% dextrose

---



### **Extraneal®**

Solution for the Long Dwell in CAPD and APD Patients, using Icodextrin instead of Glucose as the Primary Osmotic Agent.

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

---



### **Physioneal® 40**

A Solution that provides a Physiologically Balanced Environment for peritoneal membrane.

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15mmol/L, Bicarb 25mmol/L

Available in 1.36, 2.27 and 3.86% glucose

---

Dianeal® PD4, PD101 and Physioneal® 40 Prescribing Information available upon request.  
Extraneal® (7.5% Icodextrin) Product Monographs available upon request.

# Peritoneal Dialysis

## CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

### TWIN BAG®

The system of choice for Canadian CAPD patients.

When it comes to your patients, trust TwinBag

- Simple Design, Successful Therapy
- Engineered for optimal flush efficiency
- Protects against touch contamination
- Allows effective removal of potential contaminants
- Designed to minimize fibrin blockage
- Latex free



CATALOGUE NUMBER	CONCENTRATION	FILL VOLUME (mL)	DRAIN BAG VOLUME (mL)	PACK FACTOR
------------------	---------------	------------------	-----------------------	-------------

#### DIANEAL® - PD4

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

JB9415	1.50% Dextrose	1500	2000	6/case
JB9425	2.50% Dextrose	1500	2000	6/case
JB9435	4.25% Dextrose	1500	2000	6/case
JB9416	1.50% Dextrose	2000	2500	6/case
JB9426	2.50% Dextrose	2000	2500	6/case
JB9436	4.25% Dextrose	2000	2500	6/case
JB9418	1.50% Dextrose	2500	3000	4/case
JB9428	2.50% Dextrose	2500	3000	4/case
JB9438	4.25% Dextrose	2500	3000	4/case

#### DIANEAL® - PD101

Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L


JB9806	0.50% Dextrose	2000	2500	6/case
JB9816	1.50% Dextrose	2000	2500	6/case
JB9826	2.50% Dextrose	2000	2500	6/case
JB9836	4.25% Dextrose	2000	2500	6/case

# Peritoneal Dialysis

## CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

### TWIN BAG®

(Continued)



#### COLOUR CODED CAP CONFIRMS DEXTROSE CONCENTRATION

- White = 0.50% dextrose (Dianeal®)
- Yellow = 1.50% dextrose (Dianeal®) and 1.36% glucose (Physioneal®)
- Green = 2.50% dextrose (Dianeal®) and 2.27% glucose (Physioneal®)
- Orange = 4.25% dextrose (Dianeal®) and 3.86% glucose (Physioneal®)
- Purple = 7.5% Icodextrin (Extraneal®)

CATALOGUE NUMBER	CONCENTRATION	FILL VOLUME (mL)	DRAIN BAG VOLUME (mL)	PACK FACTOR
------------------	---------------	------------------	-----------------------	-------------



#### EXTRANEAL®

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

JB9912	7.5% Icodextrin Solution	2000	2500	6/case
JB9913	7.5% Icodextrin Solution	2500	3000	4/case






#### PHYSIONEAL®

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25mmol/L

JB8249	1.36% Glucose	2000	2500	6/case
JB8251	2.27% Glucose	2000	2500	6/case
JB8254	3.86% Glucose	2000	2500	6/case
JB8250	1.36% Glucose	2500	3000	5/case
JB8252	2.27% Glucose	2500	3000	5/case


# Peritoneal Dialysis

## CAPD – CONTINUOUS AMBULATORY PERITONEAL DIALYSIS

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
<b>Disconnect Y-Sets</b>			
5C4366P	CAPD Disconnect Y-Set, Bulk	60/case	
<b>Miscellaneous CAPD Supplies</b>			
5C4462	Empty Sterile Bag, 3L, Spike	48/case	
XMC4284	Empty Sterile Heater/Drain Bag, 3L, Luer Lock	30/box	

# Peritoneal Dialysis

## CAPD Solution Warmer

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
JBWarmerU	CAPD Solution Warmer	1 each	

Designed to provide patients with simple and efficient system for warming TwinBag solutions.

- Ensures integrity of PD solutions
- Preset to a comfortable temperature
- Secondary overheat thermal protector shutoff
- Lightweight and portable
- Requires only 20 watts of electricity
- Made in Canada

# Peritoneal Dialysis

## APD – AUTOMATED PERITONEAL DIALYSIS

### HomeChoice Claria

- Initial Drain (I-Drain) Bypass - Maximizing patient comfort and confidence
- Sharesource 2-way connectivity platform for Remote Patient Management (RPM)
- 2-line OLED screen that eliminates alternating messages and helps improve the user experience (1)
- Increased number of languages to 38



HomeChoice Claria

#### DEVICES

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
5C6M40	HOMECHOICE CLARIA APD	1 each
5C6M40R	HOMECHOICE CLARIA APD REFURB	1 each
5C7346	HOMECHOICE CLARIA CARRYINGCASE	1 each
5CGM01	CLARIA SHARESOURCE PLATFORM	1 each
5C6M70	4G KIT CELLULAR MODEM AND POWER	1 each
BXU005345	POWER CORD FOR CLARIA 125V NEMA	1 each
0719003471A1*	HOMECHOICE CLARIA PATIENT AT HOME GUIDE (PAHG)	1 each
0719001577C1**	HOMECHOICE CLARIA Clinician guide English	1 each



4G Modem



<sup>1)</sup> Baxter International Inc. Deerfield, IL

\*0719003471A1FRE: PAHG French. \*\*0719001577C1FRE: Clinician guide French

# Peritoneal Dialysis

## AMIA Automated Peritoneal Dialysis System with Sharesource Connectivity Platform

### AVAILABLE FOR ADULT PATIENTS

- 40% smaller and 30% lighter than the HOMECHOICE PRO Automated PD System
- Provides a user-friendly PD experience that guides patients through a treatment or therapy
- Equipped with voice guidance, touch screen, animations, and remote management capabilities



### DEVICES

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
5C9320	AMIA APD SYSTEM V2.0	1
5C9320R	AMIA APD SYSTEM V2.0 Refurb	1
5C6012	Country Kit, V 2.0, Canada, English	1
5C9346	Cycler Travel Case	1
5C9400	AMIA Sharesource Platform	NA
R5C5522	Personal Data Key 5-Pack	5
5C6M70**	4G Kit cellular Modem and Power	1
5C5479	Set 4 prong. APD Cycler	30

\*\*5C6M70: distribution started Feb 2021. It is replacing 5C6M20



# Peritoneal Dialysis

## APD – AUTOMATED PERITONEAL DIALYSIS

### HomeChoice Pro™ and HomeChoice™ Automated Peritoneal Dialysis Systems

**AVAILABLE FOR  
PEDIATRIC TO ADULT  
PATIENTS**

- Pneumatic pumps – no scales!
- No re-calibration necessary after moving the machine
- Patient management with ProCard and RenalSoft PD

The low-fill mode option for pediatric and low-fill volume adult patients offers:

- Fill volume range from 60 – 3000 mL
- Programmable minimum drain time and volume parameters for more closely regulated therapy
- Alarm limits and programmable positive and negative UF alarms ensure proper drainage and increase treatment flexibility





**HomeChoice Pro™**

#### DEVICES

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR
5C4471RU*	HomeChoice™ Automated PD System	1 each
5C4471RUL	HomeChoice™ Automated PD System with Mandarin Language	1 each
5C8310RU**	HomeChoice PRO™ Automated PD System (110 volt)	1 each

## Peritoneal Dialysis Software

- APD therapies are written by the clinician and transferred to the cyclor with the PRO Card
- Window's based software enables increased monitoring of patient adherence, drain function and ultrafiltration

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
R5C4522	HomeChoice PRO™ Card	6/case	
4994000449	Homechoice Pro USB SmartMedia Card Reader	1 each	

\*5C4471RU is replacing 5C4471R. \*\*5C8310RU is replacing 5C8310P



# Peritoneal Dialysis

## APD – AUTOMATED PERITONEAL DIALYSIS

### APD Luer Lock Solutions

CATALOGUE NUMBER	CONCENTRATION	FILL VOLUME (mL)	PACK FACTOR
------------------	---------------	------------------	-------------

#### DIANEAL® - PD4

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 40 mmol/L

JB4769L	1.50% Dextrose	3000	4/case
JB4770L	2.50% Dextrose	3000	4/case
JB4771L	4.25% Dextrose	3000	4/case
JB9719L	1.50% Dextrose	5000	2/case
JB9729L	2.50% Dextrose	5000	2/case
JB9749L	4.25% Dextrose	5000	2/case

#### DIANEAL® - PD101

Na 132, Ca 1.62, Mg 0.75, Cl 101.75, Lactate 35 mmol/L

JB4821L	0.50% Dextrose	2000	6/case
JB5618L	1.50% Dextrose	2500	5/case
JB5658L	2.50% Dextrose	2500	5/case
JB5619L	1.50% Dextrose	5000	2/case
JB5659L	2.50% Dextrose	5000	2/case
JB5669L	4.25% Dextrose	5000	2/case

#### EXTRANEAL®

Na 132, Ca 1.75, Mg 0.25, Cl 96, Lactate 40 mmol/L

JB9923LP	7.5% Icodextrin Solution	2500	5/case
----------	--------------------------	------	--------

#### PHYSIONEAL® 40

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25 mmol/L

JB8291	1.36% Glucose	2500	5/case
JB8290	2.27% Glucose	2500	5/case
JB8289	3.86% Glucose	2500	5/case






#### PHYSIONEAL® 40

Na 132, Ca 1.25, Mg 0.25, Cl 95, Lactate 15, Bicarb 25 mmol/L


CJPE8280	1.36% Glucose	5000	2/case
CJPE8281	2.27% Glucose	5000	2/case
CJPE8282	3.86% Glucose	5000	2/case

# Peritoneal Dialysis





## APD – AUTOMATED PERITONEAL DIALYSIS

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
<b>APD DISPOSABLE LUER SETS</b>			
R5C4479C	4-Prong Cyclor Set with Cassette	30/case	
R5C4478	8-Prong Cyclor Set with Cassette	30/case	
R5C4427	5-Prong Manifold Set	30/case	
N5C8305C	Low Recirculation Volume APD Set  Low Recirculation Volume APD Set has in-line medication port, Luer Lock to bag connectors, a shrouded patient connector which covers all transfer set threads, and a smaller internal diameter and shorter patient line length which reduces patient line recirculation volume to 17mL for a more efficient therapy.	30/case	
<b>APD EXTENSION LINE</b>			
5C4480C	12 Foot Patient Extension Line, Easy Lock™ Connector (3.6M)	30/case	

# Peritoneal Dialysis

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
<b>APD DRAIN SUPPLIES</b>			
XMC4284	Empty Sterile Heater\Drain Bag, 3L, Luer Lock	30/box	
5C4462	Empty Sterile Bag, 3L, Spike	48/case	
5C4476	Effluent Sample Bag, 150ml, Luer	6/case	
5C4145P	Cycler Drainage Bag Set (2 x 15L per pack), Spike	15 pk/case	
5C4512	APD Drain Manifold (Connects 2 15 Litre drain bags), Spike	30/case	
5C4464P	12 Foot Drain Line Extension Set (3.6M), Spike	30/case	






# Peritoneal Dialysis

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
<b>APD CARRYING CASES</b>			
5C7346	Homechoice Claria Carrying Case	1 each	
5C9346	AMIA Cyclor Travel Case	1 each	
5C4406	HomeChoice: Soft-sided Carrying Case	1 each	
5C8346	HomeChoice: Shipping Case with wheels	1 each	
<b>CARTS</b>			
EK8029	Home Choice Cart with sliding shelf & stationary arms Wood melamine white cart with 4 plastic casters, hook, 1 shelf, 1 drawer. 3 side raise panel made in maple wood These unit are pack individually in boxes with Styrofoam protective padding	1 each	

# Peritoneal Dialysis

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
<b>Transfer Sets</b>			
<b>-Disconnect System</b>			
5C4483	MiniCap Extended Life Transfer Set with Twist Clamp (10cm)*	6/case	
<b>5C4482</b>	<b>MINICAP EXTENDED LIFE TRANSFER SET WITH TWIST CLAMP (15CM)*</b>	<b>6/case</b>	
R5C4484	MiniCap Extended Life Transfer Set with Twist Clamp (22cm)*	6/case	
5C4449	MiniCap Extended Life Transfer Set with Roller Clamp (15cm)*	6/case	
<b>-Non-Disconnect System</b>			
<b>Transfer Set Adapters</b>			
5C4136	Luer-Lock to Spike Adapter Set	12/case	
<b>Catheter Supplies</b>			
5C4129	Locking Titanium Adapter (2 pieces)	5/bx,24bx/ca	
5C4169	Locking Cap for 5C4129 Titanium Adapter	25/case	
PC4171	Catheter Clamp	12/case	
<b>5C4319</b>	<b>IMMOBILE™ CATHETER AND TRANSFER SET STABILIZING DEVICE</b>	<b>30/case</b>	
5C4320	IMMOBILE A/C, CATHETER STABILIZER	30/case	

# Peritoneal Dialysis

CATALOGUE NUMBER	DESCRIPTION	PACK FACTOR	PICTURE
<b>Disconnect Caps with Pre-filled Povidone-Iodine Solution</b>			
<b>5C4466P</b>	<b>MiniCap Disconnect Cap</b> (English labeling only)	60/case	
SPC4466	MiniCap Disconnect Cap (same as RPC4466)	60/case	
SPC4456	Flexicap Disconnect Cap	30/case	
<b>Clamps</b>			
5C4527	Outlet Port Clamp, Short-Nosed, Blue	12/case	
5C4957	UltraClamp, Red (Not for use as outlet port clamp)	12/case	
<b>EZ-Aide™</b>			
5C4505	Assist Device for CAPD or APD patients	1 each	

**\*Bolded codes are preferred products**